

UNIVERSITY OF LEIPZIG

ADVANCED LABS

## Lab report

# Doppler-free Rb saturation spectroscopy with an external cavity diode laser

Jamal Ghaith 3792970

Anas Roumieh 3766647

Conducted on: 14.05.2024

# Contents

<b>1</b>	<b>Introduction</b>	<b>1</b>
<b>2</b>	<b>Analysis</b>	<b>2</b>
2.1	Task 1 . . . . .	2
2.1.1	Scaling the data . . . . .	2
2.1.2	Determining the finesse . . . . .	3
2.2	Task 2 . . . . .	4
2.3	Task 3 . . . . .	6
2.4	Task 4 . . . . .	7
<b>3</b>	<b>Conclusion</b>	<b>8</b>
	<b>Appendices</b>	<b>9</b>
	<b>Bibliography</b>	<b>10</b>

# 1 Introduction

$$\mathcal{F} = \frac{\nu_{\text{FSR}}}{\delta\nu} \tag{1}$$

## 2 Analysis

### 2.1 Task 1

We were instructed to scale our measurement data using the FPI peaks in addition to determining the finesse.

#### 2.1.1 Scaling the data

From [1], we know that our FSR is 1 GHz. Therefore, if the average spacing between peaks is calculated, we can determine the conversion factor and scale our data accordingly.




Figure 1: FPI peaks. 1a Raw data. 1b Flipped and normalized data, with peaks highlighted

The average spacing between peaks was calculated to be  $\approx 56.15$ , meaning there are

$$\approx \frac{1 \text{ GHz}}{56.15 \text{ counts}} \approx 0.0178 \frac{\text{GHz}}{\text{count}}$$

Scaling the data using this conversion factor, we obtain the following plot:




Figure 2: Scaled data using the FPI peaks.

### 2.1.2 Determining the finesse

Using equation 1, the FWHM for a selected FPI peak can be used to find the finesse. The following shows a Lorentzian fit on a selected peak:




Figure 3: Lorentzian fit on a selected FPI peak.

From the fit,  $\delta\nu \approx 0.00476$  GHz Hence, the finesse is

$$\mathcal{F} = \frac{\nu_{FSR}}{\delta\nu} = \frac{1 \text{ GHz}}{0.00476 \text{ GHz}} \approx 209.9$$

## 2.2 Task 2

First, we assign the peaks in the rubidium spectrum to the correct isotopes in specific ground states. The gas cell contains a mixture of  $^{85}\text{Rb}$  and  $^{87}\text{Rb}$ . The  $D_2$  line (780.2 nm) in  $^{85}\text{Rb}$  corresponds to the transition from the excited state  $5^2P_{3/2}$  (with  $F = 1, 2, 3, 4$ ) to the ground state  $5^2S_{1/2}$  (with  $F = 2, 3$  possible). For  $^{87}\text{Rb}$  the  $D_2$  line corresponds to the same transition but with  $F = 0, 1, 2, 3$  possible in the excited state and  $F = 1, 2$  in the ground state.

Now, we want to analyze the effect of changing the temperature and injection current of the diode. Hence we show plots of the saturation spectra for constant temperature and varied current and vice versa. We discuss in the following the behaviour observed.

The variation of temperature and current could lead to the following effects:

1. Mode hopping:

Mode hopping occurs in lasers when the laser's output frequency abruptly shifts between different longitudinal modes of the laser cavity. Which could happen due to changes in conditions such as temperature and power (injection current) which alter the wavelength of maximal gain. When mode hopping happens, a mixture of different frequencies might be emitted before a new mode completely dominates the optical power of the diode. This can result in multiple distinct peaks and patterns appearing at various temperatures and injection currents. Mode hopping leads to a change in the frequency of the beam and hence the intensity of the spectrum.

2. Size change of the active region in the laser diode:

An increase in the temperature or injection current (which also increases temperature) will cause the semiconductor material which make up the laser to expand. This affects the dimensions of the active region which changes the optical modes in laser cavity.

3. Spectral broadening:

First, we discuss the constant temperature case. Sub-figure (f) shows a significant change upon changing the current from  $I = 121\text{ mA}$  to  $I = 122\text{ mA}$  with a drop in the intensity, which could be interpreted as a mode jump.




Figure 4: Saturation spectra for variation of injection currents

For the case of constant injection current, a mode jump occurs between  $T = 19.9^\circ\text{C}$  and  $T = 20.1^\circ\text{C}$ . We notice that changing the current affects the spectrum more drastically. This could be due to the fact that a slight change in the current leads to more heating compared to the temperature increments we used.




Figure 5: Saturation spectra for variation of temperature

**2.3 Task 3**



**2.4 Task 4**

### 3 Conclusion

# Appendices

## Bibliography

- [1] F. Jung, “Doppler-free rb saturation spectroscopy,” 2018.